Biotecnologia

ACH5545 Engenharia Genética

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2º Semestre 2024

Docente:

Felipe Chambergo (fscha@usp.br) - <https://sites.usp.br/lbbp>

Monitores:

Augusto Roldan Gonçalves - augusto.roldan@usp.br

Henrique dos Santos Hernandes - hernandesrique@usp.br

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Tec. Pedro Manoel dos Santos - pedroms@usp.br

Créditos: 4

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Bibliografia básica

1. Ausubel, Frederick M, Roger Brent, Robert E. Kingston, David D. Moore, J.G. Seidman, John A. Smith. Current Protocols in Molecular Biology. 2003. Kevin Struhl (eds.) John Wiley & Sons, Inc.
 2. Dixit, R; Bisen, K; Kumar, A; Borah, A e Keswani, C. .LAB MANUAL ON MOLECULAR BIOLOGY, First Edition : 2016. ISBN : 978-81-909182-7-5, Published by : Media Associates J-281/86, Kartar Nagar, Delhi-110053, India.
 3. Green, Michael R., Howard Hughes Medical Institute, University of Massachusetts Medical School; Joseph Sambrook, Peter MacCallum Cancer Institute, Melbourne, Australia. Molecular Cloning: A Laboratory Manual (Fourth Edition). 2012 • 2,028 pp., ISBN 978-1-936113-42-2
 4. Wilson, K and John Walker (Eds). Principles and techniques of biochemistry and molecular biology. 2010. 7th ed. Cambridge University Press The Edinburgh Building, Cambridge CB2 8RU, UK. ISBN 978-0-521-51635-8 (hardback)
 5. Molecular Cloning, A Laboratory Manual, 4th Edition, www.molecularcloning.org. Cold Spring Harbor Protocols, www.cshprotocols.org (WEBSITE).
- Leitura de artigos.
 - Protocolos fornecidos pelos fabricantes dos kits e reagentes utilizados nas atividades práticas.

MOLECULAR CLONING NEW EDITION

A Laboratory Manual


By Michael R. Green and Joseph Sambrook

Cold Spring Harbor Laboratory Press

A Book Preview Site for Molecular Cloning 4th Edition

Molecular Cloning: A Laboratory Manual (Fourth Edition)

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






SnapGene Viewer

<https://www.snapgene.com/>



Beyond the Basics of Molecular Cloning

SnapGene is the most popular cloning tool for a reason. It's fast, smart and extremely user-friendly

- Intuitive technology identifies design flaws in cloning procedures so they can be corrected
- Simulate standard PCR using your own primers, or allow SnapGene to design them automatically
- Specialised cloning tools ensure fast accurate construct design for all major molecular cloning techniques

GenBank: MT126808 SARS-CoV-2 - Brasil

SS - signal sequence
 NTD - N-terminal domain
 RBD - receptor-binding domain
 SD1 - subdomain 1
 SD2 - subdomain 2
 S1/S2 = S1/S2 protease cleavage site
 S2' = S2' protease cleavage site

FP = fusion peptide
 HR1 = heptad repeat 1
 CH = central helix
 CD = connector domain
 HR2 = heptad repeat 2
 TM = transmembrane domain
 CT = cytoplasmic tail

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


Open Reading Frames (Orf)

Non-Structural Proteins (Nsp)

| Experimentally-determined SARS-CoV-2 structure with PDB ACCESS CODE provided
| De Novo model
| Homology model
| © Denotes N7-methyl Guanine Cap
| AAAAA Denotes Poly-A Tail

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Protein	Mol. weight (kDa)	Seq. similarity with SARS-CoV	Description
Nsp1	19.8	91.1%	Suppresses host antiviral response
Nsp2	70.5	82.9%	
Nsp3	217.3	86.5%	Nsp3-Nsp4-Nsp6 complex involved in viral replication
Nsp4	56.2	90.8%	Nsp3-Nsp4-Nsp6 complex involved in viral replication
Nsp5	33.8	98.7%	Main protease (3C-like)
Nsp6	33.0	94.8%	Nsp3-Nsp4-Nsp6 complex involved in viral replication
Nsp7	9.2	100.0%	Nsp7-Nsp8 complex is part of RNA polymerase
Nsp8	21.9	99.0%	Nsp7-Nsp8 complex is part of RNA polymerase
Nsp9	12.4	98.2%	ssRNA binding
Nsp10	14.8	99.3%	Essential for Nsp16 methyltransferase activity
Nsp11	1.3	92.3%	Short peptide
Nsp12	106.7	98.3%	RNA polymerase
Nsp13	66.9	100.0%	Helicase/triphosphatase
Nsp14	59.8	98.7%	3'-5' exonuclease
Nsp15	38.8	95.7%	Uridine-specific endoribonuclease
Nsp16	33.3	98.0%	RNA-cap methyltransferase
S	141.2	87.0%	Spike protein, mediates binding to ACE2
Orf3a	31.1	85.1%	Activates the NLRP3 inflammasome
Orf3b	6.5	9.5%	
E	8.4	96.1%	Envelope protein, involved in virus morphogenesis and assembly
M	25.1	96.4%	Membrane glycoprotein, predominant component of the envelope
Orf6	7.3	85.7%	Type I IFN antagonist
Orf7a	13.7	90.2%	Virus-induced apoptosis
Orf7b	5.2	84.1%	
Orf8	13.8	45.3%	
N	45.6	94.3%	Nucleocapsid phosphoprotein, binds to RNA genome
Orf9b	10.8	84.7%	Type I IFN antagonist
Orf9c	8.0	78.1%	
Orf10	4.4	-	

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
Proteína S

Schematic of 2019-nCoV S primary structure colored by domain. Domains that were excluded from the ectodomain expression construct or could not be visualized in the final map are colored white. SS, signal sequence; S2', S2' protease cleavage site; FP, fusion peptide; HR1, heptad repeat 1; CH, central helix; CD, connector domain; HR2, heptad repeat 2; TM, transmembrane domain; CT, cytoplasmic tail. Arrows denote protease cleavage sites.

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