Biotecnologia

ACH5545 Engenharia Genética

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2º Semestre 2024

Docente:

Felipe Chambergo (fscha@usp.br) - <https://sites.usp.br/lbbp>

Monitores:

Augusto Roldan Gonçalves - augusto.roldan@usp.br

Henrique dos Santos Hernandes - hernandesrique@usp.br

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Tec. Pedro Manoel dos Santos - pedroms@usp.br

Créditos: 4

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[User Programs](#)

[News & Publications](#)

A BERKELEY LAB NATIONAL USER FACILITY


VEGA Viral EcoGenomics & Applications

April 8 - May 13, 2021

Shannon Bennett
Aude Bernheim
Edward Holmes
Mart Krupovic

Jens Kuhn
Karen Maxwell
Monir Moniruzzaman

David Pride
Simon Roux
Rachel Whitaker
Natalya Yutin

Don't miss our VEGA Symposium!

Latest News

FEBRUARY 22, 2021

[Green Algae Reveal One mRNA Encodes Many Proteins](#)

FEBRUARY 8, 2021

[An Age of CRAGE: Advances in Rapidly Engineering Non-model Bacteria](#)

JANUARY 27, 2021

[Fields of Breeders' Dreams: A Team Effort Toward Targeted Crop Improvements](#)

[LATEST JGI NEWS RELEASES](#)

Announcements

Read our [Statement on the Use of Genomics Data](#).

Our Science

How we're contributing to bioenergy research.

[LEARN MORE](#)

User Programs

Everything you need to know about getting a project started.

[LEARN MORE](#)

<https://jgi.doe.gov/>

NCBI Resources How To Sign in to NCBI

NCBI National Center for Biotechnology Information

COVID-19 Public health information

Notice of S April 24, 2022 the eRA website

Ending Structural Racism nih.gov/ending-structural-racism

All Databases

- Gene
- Genome
- GEO DataSets
- GEO Profiles
- GTR
- HomoloGene
- Identical Protein Groups
- MedGen
- MeSH
- NCBI Web Site
- NLM Catalog
- Nucleotide
- OMIM
- PMC
- PopSet
- Protein
- Protein Clusters
- Protein Family Models
- PubChem BioAssay
- PubChem Compound

Information (NIH) | SARS-CoV-2 data (NCBI) | Prevention and treatment information (HHS) | Español

Please note that eRA systems will be undergoing scheduled maintenance from 10am until 10pm Eastern US time on dependent services such as MyBibliography and Grant Reporting may be unavailable. More information is available on

and structural racism and achieve racial equity in the biomedical research enterprise.

NCBI National Center for Biotechnology Information advances science and health by providing access to genomic information.

About the NCBI | Mission | Organization | NCBI News & Blog

Submit
Deposit data or manuscripts into NCBI databases

Download
Transfer NCBI data to your computer

Learn
Find help documents, attend a class or watch a tutorial

Develop
Use NCBI APIs and code libraries to build applications

Analyze
Identify an NCBI tool for your data analysis task

Research
Explore NCBI research and collaborative projects

Popular Resources

- PubMed
- Bookshelf
- PubMed Central
- BLAST
- Nucleotide
- Genome
- SNP
- Gene
- Protein
- PubChem

NCBI News & Blog

Data for SARS-CoV-2 variants now available at NCBI

23 Apr 2021
Looking for genomes for the B.1.1.7


NCBI Home
Resource List (A-Z)
All Resources

- Chemicals & Bioassays
- Data & Software
- DNA & RNA
- Domains & Structures
- Genes & Expression
- Genetics & Medicine
- Genomes & Maps
- Homology
- Literature
- Proteins
- Sequence Analysis
- Taxonomy
- Training & Tutorials

In Focus

Comment on recent variants and spike protein changes

As seen on many occasions before, mutations are naturally expected for viruses and are most often simply neutral regional markers useful for contact tracing. The changes seen have rarely affected viral fitness and almost never affected clinical outcome but the detailed effects of these mutations remain to be determined fully. Changes in the spike protein have relevance for potential effects on both host receptor as well as antibody binding with possible consequences for infectivity, transmission potential and antibody and vaccine escape. Actual effects need to be measured and verified experimentally.

As has become evident, these few S gene mutations and some deletions are found in multiple genomic contexts (different clades in different countries) that may be an early indication for some potential advantage for these viruses but needs to be verified and does not necessarily mean change in clinical severity or transmission efficiency. [> read more](#)

Genomic epidemiology of hCoV-19

hCoV-19 data sharing via GISAID

1,225,793
submissions


Enabled by data shared via GISAID

SEARCH AND SIMILARITY QUERY

Status of Detection Systems

University of Turin (Italy)
Diagnostic detection of 2019-nCoV by real-time

hCoV-19 Submission Tracking

hCoV-19 Tracking of Variants

Public-Private Partnerships of the GISAID Initiative

The GISAID Initiative involves public-private-partnerships between the Initiative's administrative arm Freunde of GISAID e.V., a registered non-profit association, and governments of the [Federal Republic of Germany](#), the official host of the GISAID platform, [Singapore](#) and the [United States of America](#), with support from private and corporate philanthropy.

hosted by the
Federal Republic
of Germany

GISAID Resources

Free Access Credentials
Register here and join thousands of researchers around the globe.

Frequently Asked Questions
Answers about the GISAID Initiative

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ExPASy 🇨🇭 Translate

Programmatic access ↓

Translate is a tool which allows the translation of a nucleotide (DNA/RNA) sequence to a protein sequence.

DNA or RNA sequence

Please enter a DNA or RNA sequence - numbers and blanks are ignored

Output format

- Verbose: Met, Stop, spaces between residues
- Compact: M, -, no spaces
- Includes nucleotide sequence
- Includes nucleotide sequence, no spaces

DNA strands

forward reverse

Genetic codes - See NCBI's genetic codes

Standard

reset TRANSLATE!

 ExPASy is operated by the SIB Swiss Institute of Bioinformatics | Terms of Use [Back to the top](#)

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Ferramenta: BLAST - NCBI

The screenshot shows the NCBI BLAST website. At the top left is the BLAST logo, and at the top right are navigation links: Home, Recent Results, Saved Strategies, and Help. The main content area is divided into several sections:

- Basic Local Alignment Search Tool:** A descriptive paragraph explaining that BLAST finds regions of similarity between biological sequences by comparing nucleotide or protein sequences to sequence databases and calculating statistical significance. A "Learn more" link is provided.
- NEWS:** A news box with a vertical "NEWS" label on the left. The text reads: "A new feature was added to Primer-BLAST. We now offer the ability for user to run primer-blast from NCBI assembly page..". The date and time are "Tue, 23 Feb 2021 12:00:00 EST", and there is a link for "More BLAST news...".
- Web BLAST:** A section with three main options:
 - Nucleotide BLAST:** Represented by a green box with a DNA double helix icon. The text below the icon says "nucleotide ► nucleotide".
 - blastx:** Represented by a blue arrow pointing right. The text inside says "translated nucleotide ► protein".
 - tblastn:** Represented by a blue arrow pointing left. The text inside says "protein ► translated nucleotide".
 - Protein BLAST:** Represented by a blue box with a protein ribbon icon. The text below the icon says "protein ► protein".
- BLAST Genomes:** A section with a search input field containing the placeholder text "Enter organism common name, scientific name, or tax id" and a "Search" button. Below the input field are four buttons: "Human", "Mouse", "Rat", and "Microbes".

[← Edit Search](#)[Save Search](#)[Search Summary ▾](#)[? How to read this report?](#)[▶ BLAST Help Videos](#)[↶ Back to Traditional Results Page](#)

Job Title	1
RID	3T2J2RAV016 <small>Search expires on 03-25 23:12 pm</small> Download All ▾
Program	BLASTX ? Citation ▾
Database	nr See details ▾
Query ID	lcl Query_770492
Description	1
Molecule type	dna
Query Length	2320
Other reports	?

Filter Results

Organism only top 20 will appear exclude[+ Add organism](#)

Percent Identity

E value

Query Coverage

 to to to [Filter](#)[Reset](#)

Descriptions

[Graphic Summary](#)[Alignments](#)[Taxonomy](#)

Sequences producing significant alignments

[Download](#) ▾[Select columns](#) ▾

Show

100 ▾

[?](#) select all 100 sequences selected[GenPept](#)[Graphics](#)

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
<input checked="" type="checkbox"/>	glycoside hydrolase family 6 [Trichoderma reesei QM6a]	Trichoderma reesei QM6a	591	740	62%	0.0	91.20%	471	XP_006962580.1
<input checked="" type="checkbox"/>	Chain A_CELLOBIOHYDROLASE II CORE PROTEIN [Trichoderma reesei]	Trichoderma reesei	591	677	51%	0.0	91.20%	365	3CBH_A
<input checked="" type="checkbox"/>	Chain A_CELLOBIOHYDROLASE CEL6A (FORMERLY CALLED CBH II) [Trichoderma reesei]	Trichoderma reesei	591	674	51%	0.0	91.20%	363	1QK0_A
<input checked="" type="checkbox"/>	cellobiohydrolase II [synthetic construct]	synthetic construct	590	738	61%	0.0	91.20%	447	AER26911.1
<input checked="" type="checkbox"/>	cellbiohydrolase II [Trichoderma reesei]	Trichoderma reesei	590	739	62%	0.0	90.91%	471	ADC83999.1
<input checked="" type="checkbox"/>	Chain A_CELLOBIOHYDROLASE II [Trichoderma reesei]	Trichoderma reesei	589	675	51%	0.0	90.91%	365	1CB2_A
<input checked="" type="checkbox"/>	hypothetical protein TgHK011_002311 [Trichoderma gracile]	Trichoderma gracile	589	773	65%	0.0	91.20%	471	KAH0490860.1

glycoside hydrolase family 6 [Trichoderma reesei QM6a]

NCBI Reference Sequence: XP_006962580.1

[Download Protein](#) [FASTA](#) [Align/BLAST](#)[Go to:](#) [⌕](#)

LOCUS XP_006962580 471 aa linear PLN 05-FEB-2020
 DEFINITION glycoside hydrolase family 6 [Trichoderma reesei QM6a].
 ACCESSION XP_006962580
 VERSION XP_006962580.1
 DBLINK BioProject: [PRJNA225530](#)
 BioSample: [SAMN02746107](#)
 DBSOURCE REFSEQ: accession [XM_006962518.1](#)
 KEYWORDS RefSeq.
 SOURCE Trichoderma reesei QM6a

ORGANISM *Trichoderma reesei QM6a*
 Eukaryota; Fungi; Dikarya; Ascomycota; Pezizomycotina;
 Sordariomycetes; Hypocreomycetidae; Hypocreales; Hypocreaceae;
 Trichoderma.
 REFERENCE 1 (residues 1 to 471)
 AUTHORS Martinez, D., Berka, R.M., Henriessat, B., Saloheimo, M., Arvas, M., Baker, S.E., Chapman, V., Chernikov, G., Coutinho, P.M., Cullen, D., Danchin, E.S., Grigoriev, I.V., Harris, F., Jackson, M., Kubicek, C.P., Han, C.S., Ho, I., Laxroo, L.F., de Leon, A.L., Magnuson, J.K., Merino, S., Mira, M., Nelson, B., Putnan, N., Robertse, B., Salamov, A.A., Schell, M., Terry, A., Trayer, M., Westerholm-Barvinen, A., Schoch, C.L., Yao, J., Barabote, R., Nelson, M.A., Detter, C., Bruce, D., Kuske, C.R., Xie, G., Richardson, F., Rokhsar, D.S., Lucas, S.M., Rubin, E.M., Dunn-Coleman, N., Ward, M. and Brettin, T.S.
 TITLE Genome sequencing and analysis of the biomass-degrading fungus *Trichoderma reesei* (syn. *Hypocrea jecorina*)
 JOURNAL Nat. Biotechnol. 26 (5), 553-560 (2008)
 PUBMED [18454138](#)
 REMARK Erratum [Nat Biotechnol. 2008 Oct;26(10):1193.. Barbote, Ravi [corrected to Barabote, Ravi]]
 REFERENCE 2 (residues 1 to 471)
 CONSTRM NCBI Genome Project
 TITLE Direct Submission
 JOURNAL Submitted (01-FEB-2020) National Center for Biotechnology

FEATURES
 Location/Qualifiers
 source 1..471
 /organism="Trichoderma reesei QM6a"
 /strain="QM6a"
 /type_material="type material of Trichoderma reesei"
 /db_xref="taxon:431251"
 /chromosome="Unknown"
 Protein 1..471
 /product="glycoside hydrolase family 6"
 /calculated_mol_wt=49522
 Region 30..62
 /region_name="fCBD"
 /note="Fungal-type cellulose-binding domain: smart00236"
 /db_xref="CDD:197593"
 Region 132..434
 /region_name="Glyco_hydro_6"
 /note="Glycosyl hydrolases family 6; pfam01341"
 /db_xref="CDD:334499"
 CDS 1..471
 /locus_tag="TRIREDRAFT_72567"
 /coded_by="XM_006962518.1:53..1468"
 /note="cellulohydrolase II (Cel16A); N-terminal CBM1 module"
 /db_xref="GeneID:18488147"
 /db_xref="InterPro:IPR000254"
 /db_xref="InterPro:IPR001524"
 /db_xref="InterPro:IPR001524"

ORIGIN
 1 mivgilttla tlaclaasvp leergacsav wqggggqns gptccasgt cyvandyayq
 61 clpqaasasa straaatstr vsptrssaa atpppggett rrvpvgagta tyeqgnfvqv
 121 tpenayyas evsalaipal tgamataaaa vakvpsfmvl dtdktpime qladlirtan
 181 knngnyaqgf vvyldprdc aalaangeys iadggvakyk nydrtirqlv veyadrtll
 241 vlepdelanl vmlgtpkca naqsayleci nyavtqinlp nvamylagah agwlgvpaq
 301 dpaqlfanv ykmaaspral rqlatnvany ngwnitppp ytgqnavyne klyibaqlpl
 361 lanhgwanaf fitdqgrgk qptggqvgd wenvigtqfg lrpsantqds lldefvwwk
 421 ggeodgteds saprfdhca lpdalqpapq agawfgyfv qlitnanpef 1
 //

Change region shown ▾

Customize view ▾

Analyze this sequence ▾

Run BLAST

Identify Conserved Domains

Highlight Sequence Features

Find in this Sequence

Show in Genome Data Viewer

Articles about the TRIREDRAFT_72567 gene ▾

Domain architecture divergence leads to functional divergence in bindi [J Biol Chem. 2020]

Genome sequencing and analysis of the biomass-degrading fungus [Nat Biotechnol. 2008]

See all...

Pathways for the TRIREDRAFT_72567 gene ▾

Metabolic pathways

Starch and sucrose metabolism

Reference sequence information ▾

RefSeq mRNA

See reference mRNA sequence for the TRIREDRAFT_72567 gene (XM_006962518.1).

PubMed (Weighted)

Related Structures (Summary)

Recent activity ▾

Turn Off Clear

glycoside hydrolase family 6 [Trichoderma reesei QM6a] Protein

Trichoderma reesei QM6a glycoside hydrolase family 6 (TRIREDRAFT_7) nucleotide

glycoside hydrolase family 7 [Trichoderma reesei QM6a] Protein

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
Ferramenta: GeneDoc

GeneDoc - [Seqs]

File Project Edit Arrange Shade Groups Score Tree Reports Plot Window Help

C S G [Icons]

C Q P E S H I L D [Icons] M U

```
gDNA : * 120 * 140 * 160 * 180 * 200
       CCGAATACCACATCGTGTGATGCTTCAGCGGCACGGAAGTACACAGACACCGCCTGTATAAAAGGGGACTGTGACCCTGTATGAGCGGCAACATGGTCT : 200
cDNA : ----- : -

gDNA : * 220 * 240 * 260 * 280 * 300
       CACAGCAGCTCACTGAAGAGGCTTGAAGATCACCTCTGTGTATTGCACCATGATTTGTGGCATTCTCACCACGCTGGCTACGCTGGCCACACTCGCA : 300
cDNA : CACAGCAGCTCACTGAAGAGGCTTGAAGATCACCTCTGTGTATTGCACCATGATTTGTGGCATTCTCACCACGCTGGCTACGCTGGCCACACTCGCA : 100

gDNA : * 320 * 340 * 360 * 380 * 400
       GCTAGTGTGCCTCTAGAGGAGCGGCAAGCTTGCTCAAGCCTCTGGTAATTATGTGAACCTCTCAAGAGACCCAAATACTGAGATATGTCAAGGGGCCAA : 400
cDNA : GCTAGTGTGCCTCTAGAGGAGCGGCAAGCTTGCTCAAGCCTCTGG-----BGCCAA : 151

gDNA : * 420 * 440 * 460 * 480 * 500
       TGTGGTGGCCAGAATTGGTCGGGTCGGACTTGCTGTGCTTCCGGAAGCACATGCGTCTACTCCAACGACTATTACTCCCAGTGTCTTCCCGGCGCTGCAA : 500
cDNA : TGTGGTGGCCAGAATTGGTCGGGTCGGACTTGCTGTGCTTCCGGAAGCACATGCGTCTACTCCAACGACTATTACTCCCAGTGTCTTCCCGGCGCTGCAA : 251

gDNA : * 520 * 540 * 560 * 580 * 600
       GCTCAAGCTCGTCCACGGCGCGCGCTCGACGACTTCTCGAGTATCCCCACAACATCCCGGTGGAGCTCCGCGACGCTCCACCTGGTTCTACTACTAC : 600
cDNA : GCTCAAGCTCGTCCACGGCGCGCGCTCGACGACTTCTCGAGTATCCCCACAACATCCCGGTGGAGCTCCGCGACGCTCCACCTGGTTCTACTACTAC : 351

gDNA : * 620 * 640 * 660 * 680 * 700
       CAGAGTACCTCCAGTCGGATCGGGAACCGCTACGTATTCAGGCCAACCTTTTGTGGGGTCACTCCTTGGGCCAATGCATATTACGCCCTCTGAAGTTAGC : 700
cDNA : CAGAGTACCTCCAGTCGGATCGGGAACCGCTACGTATTCAGGCCAACCTTTTGTGGGGTCACTCCTTGGGCCAATGCATATTACGCCCTCTGAAGTTAGC : 451

gDNA : * 720 * 740 * 760 * 780 * 800
       AGCCTCGCTATTCTAGCTTGACTGGAGCCATGCCACTGCTGCAGCAGCTGTGCGAAAGGTTCCCTCTTTATGTGGCTGTAGGTCCCTCCCGGAACCAA : 800
cDNA : AGCCTCGCTATTCTAGCTTGACTGGAGCCATGCCACTGCTGCAGCAGCTGTGCGAAAGGTTCCCTCTTTATGTGGCT----- : 531

gDNA : * 820 * 840 * 860 * 880 * 900
       GGCAATCTGTACTGAAGGCTCATCATTCACTGCAGAGATACTCTTGACAAGACCCCTCTCATGGAGCAAACCTTGGCCGACATCCGCACCGCCAACAAG : 900
cDNA : -----AGATACTCTTGACAAGACCCCTCTCATGGAGCAAACCTTGGCCGACATCCGCACCGCCAACAAG : 595

gDNA : * 920 * 940 * 960 * 980 * 1000
       AATGGCGGTAACATATGCCGACAGTTTGTGGTGTATGACTTGCCGGATCCGGATTGCGCTGCCCTTGCCCTCGAATGGCGGAATACTCTATGCCGATGGTG : 1000
cDNA : AATGGCGGTAACATATGCCGACAGTTTGTGGTGTATGACTTGCCGGATCCGGATTGCGCTGCCCTTGCCCTCGAATGGCGGAATACTCTATGCCGATGGTG : 695

gDNA : * 1020 * 1040 * 1060 * 1080 * 1100
       GCGTCGCCAAATATAAGAACTATATCGACACCAATTCGTCAAATTGTCGTGGAATATCCGATATCCGGACCCCTCCTGGTTATTGGTATGAGTTTAAACAC : 1100
cDNA : GCGTCGCCAAATATAAGAACTATATCGACACCAATTCGTCAAATTGTCGTGGAATATCCGATATCCGGACCCCTCCTGGTTATTG----- : 779
```

>File >Export >-All sequences – Clustal (aln)

CLUSTAL W(1.60) multiple sequence alignment

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[← Edit Search](#) [Save Search](#) [Search Summary](#) [? How to read this report?](#) [BLAST Help Videos](#) [↶ Back to Traditional Results Page](#)

Job Title	Protein Sequence
RID	89FZAEHH016 <small>Search expires on 04-26 09:03 am</small> Download All ▼
Program	BLASTP ? Citation ▼
Database	nr See details ▼
Query ID	lcl Query_34109
Description	None
Molecule type	amino acid
Query Length	471
Other reports	Distance tree of results Multiple alignment MSA viewer ?

Filter Results

Organism only top 20 will appear exclude

Type common name, binomial, taxid or group name

[+ Add organism](#)

Percent Identity to **E value** to **Query Coverage** to

[Filter](#) [Reset](#)

Descriptions | [Graphic Summary](#) | [Alignments](#) | [Taxonomy](#)

Sequences producing significant alignments [Download](#) New [Select columns](#) ▼ Show [?](#)

select all 100 sequences selected [GenPept](#) [Graphics](#) [Distance tree of results](#) [Multiple alignment](#) New [MSA Viewer](#)

	Description	Scientific Name	Max Score	Total Score	Query Cov.	E value	Per. Ident.	Acc. Len.	Accession
<input checked="" type="checkbox"/>	glycoside hydrolase family 6 [Trichoderma reesei QM6a]	Trichoderma ree...	959	959	100%	0.0	100.00%	471	XP_006962580.1
<input checked="" type="checkbox"/>	cellbiohydrolase II [Trichoderma reesei]	Trichoderma reesei	958	958	100%	0.0	99.79%	471	ADC83999.1
<input checked="" type="checkbox"/>	cellobiohydrolase II [Trichoderma reesei]	Trichoderma reesei	956	956	100%	0.0	99.79%	471	AAG39980.1
<input checked="" type="checkbox"/>	cellobiohydrolase II [Trichoderma viride]	Trichoderma viride	956	956	100%	0.0	99.79%	471	AAQ76094.1
<input checked="" type="checkbox"/>	unnamed [Trichoderma reesei]	Trichoderma reesei	954	954	100%	0.0	99.58%	471	AAA72922.1
<input checked="" type="checkbox"/>	GH6 Cellobiohydrolase CEL6A/CBH2 [Trichoderma parareesei]	Trichoderma par...	951	951	100%	0.0	98.73%	471	OTA06465.1
<input checked="" type="checkbox"/>	cellobiohydrolase II [Trichoderma longibrachiatum]	Trichoderma long...	939	939	100%	0.0	97.88%	470	ACZ34301.1

ProtParam

User-provided sequence:

Molecular weight: 49653.35

Theoretical pI: 5.11

Amino acid composition:

[CSV format](#)

Ala (A)	60	12.7%
Arg (R)	14	3.0%
Asn (N)	30	6.4%
Asp (D)	21	4.5%
Cys (C)	12	2.5%
Gln (Q)	21	4.5%
Glu (E)	10	2.1%
Gly (G)	40	8.5%
His (H)	4	0.8%
Ile (I)	17	3.6%
Leu (L)	37	7.9%
Lys (K)	10	2.1%
Met (M)	5	1.1%
Phe (F)	12	2.5%
Pro (P)	32	6.8%
Ser (S)	47	10.0%
Thr (T)	38	8.1%
Trp (W)	12	2.5%
Tyr (Y)	21	4.5%
Val (V)	28	5.9%
Pyl (O)	0	0.0%
Sec (U)	0	0.0%
(B)	0	0.0%
(Z)	0	0.0%
(X)	0	0.0%

Total number of negatively charged residues (Asp + Glu): 31

Total number of positively charged residues (Arg + Lys): 24

Atomic composition:

Carbon	C	2194
Hydrogen	H	3354
Nitrogen	N	594
Oxygen	O	691
Sulfur	S	17

Formula: C₂₁₉₄H₃₃₅₄N₅₉₄O₆₉₁S₁₇

Total number of atoms: 6850

Extinction coefficients:

Extinction coefficients are in units of M⁻¹ cm⁻¹, at 280 nm measured in water.

Ext. coefficient	98040
Abs 0.1% (=1 g/l)	1.974, assuming all pairs of Cys residues form cystines

Ext. coefficient	97290
Abs 0.1% (=1 g/l)	1.959, assuming all Cys residues are reduced

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CCCATCGTTCCTGTAA
GGGTAGCAAGGACATT

https://www.bioline.com/media/calculator/01_12.html

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


Plasmid Map	
NoeClone 3.0 Demo	Virtual cloning laboratory.
pDRAW32 1.1.129	Scientific software for the molecular biologist
SimVector 4.60 Demo	tool for drawing publication and vector catalog quality maps
ApE 2.0.49	A Plasmid Editor
BVTech Plasmid 5.1 Demo	DNA plasmid drawing software
XPlasMap 0.99	a DNA mapping program for Mac OS
CGView 1.0	Circular Genome Viewer
PlasmaDNA 1.4.2	A free, cross-platform PLAsmid MANipulation program
Plasmidomics 0.2	Plasmid Drawing Program.
SnapGene Viewer 3.1.2	Create, Browse, and Share richly Annotated DNA Sequence.
pLOT 1.0.10h	Plasmid mapping program
Online Tools	
PlasMapper 2.0	automatically generates and annotates plasmid maps using only the plasmid DNA sequence as input.
NetPlasmid	Online Plasmid map drawing
Savvy v0.1	Draw plasmid map tool.
EZ PLASMID MAP V1.9	Free online plasmid draw/plotter program.
WebDSV	Free Online DNA Sequence Editor and Plasmid Drawing program.

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GenBank: MT126808 SARS-CoV-2 - Brasil

SS - signal sequence
 NTD - N-terminal domain
 RBD - receptor-binding domain
 SD1 - subdomain 1
 SD2 - subdomain 2
 S1/S2 = S1/S2 protease cleavage site
 S2' = S2' protease cleavage site

FP = fusion peptide
 HR1 = heptad repeat 1
 CH = central helix
 CD = connector domain
 HR2 = heptad repeat 2
 TM = transmembrane domain
 CT = cytoplasmic tail

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