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REFERÊNCIAS

- Cox, J. & Mann, M. [MaxQuant enables high peptide identification rates, individualized p.p.b.-range mass accuracies and proteome-wide protein quantification](#). *Nature Biotechnology* 26, 1367–1372 (2008).
- Cox, J. et al. [Accurate Proteome-wide Label-free Quantification by Delayed Normalization and Maximal Peptide Ratio Extraction, Termed MaxLFQ](#). *Molecular & Cellular Proteomics* 13, 2513–2526 (2014).
- Tyanova, S., Temu, T. & Cox, J. [The MaxQuant computational platform for mass spectrometry-based shotgun proteomics](#). *Nature Protocols* 11, 2301–2319 (2016).
- Tyanova, S. et al. [Visualization of LC-MS/MS proteomics data in MaxQuant](#). *PROTEOMICS* 15, 1453–1456 (2015).
- Schaab, C., Geiger, T., Stoeck, G., Cox, J. & Mann, M. [Analysis of High Accuracy, Quantitative Proteomics Data in the MaxQB Database](#). *Molecular & Cellular Proteomics* 11, M111.014068 (2012).
- Sinitcyn, P. et al. [MaxQuant goes Linux](#). *Nature Methods* 15, 401–401 (2018).
- Tutoriais online sobre o MaxQuant oferecidos nas várias edições da MaxQuant Summer School. Disponíveis em:
<<https://www.youtube.com/channel/UCKYzYTm1cnmc0CFAMhxDO8w/videos>>

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Load Remove Write template Set experiment No fractions Set PTM

Load folder Change folder Read from file Set fractions Set parameter group Set reference channels

Experimental design file Edit experimental design

	File	Exists	Size	Data format	Parameter group	Experiment	Fraction	PTM	Reference channels
1	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO1-01.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO1-01		False	
2	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO1-02.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO1-02		False	
3	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO2-01.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO2-01		False	
4	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO2-02.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO2-02		False	
5	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO3-01.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO3-01		False	
6	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\Cold-Liver_LDLRKO3-02.raw	True	1 GB	Thermo raw...	Group 0	Cold-Liver_LDLRKO3-02		False	
7	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO1-01.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO1-01		False	
8	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO1-02.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO1-02		False	
9	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO2-01.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO2-01		False	
10	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO2-02.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO2-02		False	
11	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO3-01.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO3-01		False	
12	D:\DATA\MARIANA\240814_Testes curso analise dados BrProt\Dados sobre mouse liver RT e cold\RT-Liver_LDLRKO3-02.raw	True	1 GB	Thermo raw...	Group 0	RT-Liver_LDLRKO3-02		False	

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Group 0 Type Digestion Cross links Instrument First search

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Sequences Protein quantification Tables MS/MS analyzer Advanced

Identification Label free quantification Folder locations MS/MS fragmentation

Parameter section

Fasta files

Add Remove Change folder Identifier rule Description rule Taxonomy rule Taxonomy ID

Variation rule Test

	Fasta file path	Exists	Identifier rule	Description rule	Taxonomy rule	Taxonomy ID	Organism	Variation rule
1	D:\DATABASE\240814_uniprofkb_mus_musc...	True	%T%T%T%	%T%		10090	Mus musculus	%T%T%T%T%T%

1 item 1 selected

Include contaminants ☒

Min. peptide length 7

Max. peptide mass [Da] 6000

Min. peptide length for unspecific search 8

Max. peptide length for unspecific search 25

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